

Human VEGFR2 protein (Recombinant) (STJP000525)

STJP000525

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human VEGFR2-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant human Vascular endothelial growth factor receptor 2 (VEGFR2) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant VEGFR2 (FLK-1) can be stored in working aliquots at 2-8°C for one month, or at-20°C to-70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	3791
Gene Symbol	KDR
Uniprot ID	VGFR2_HUMAN
Immunogen Sequence	