

## Human VEGFR2 protein (Recombinant) (Fc-Tag) (STJP000526)

STJP000526

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human VEGFR2-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant KDR was lyophilized from 0.2 µm filtered PBS solution pH7.4
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage Instruction</b>	The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">3791</a>
<b>Gene Symbol</b>	<a href="#">KDR</a>
<b>Uniprot ID</b>	<a href="#">VGFR2_HUMAN</a>
<b>Immunogen</b>	
<b>Sequence</b>	