

Mouse 4-1BB Ligand (TNFSF9) protein (Recombinant) (STJP000540)

STJP000540

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse 4-1BB Ligand (TNFSF9)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse 4-1BB Ligand is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Mouse 4-1BBL can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	21950
Gene Symbol	Tnfsf9
Uniprot ID	TNFL9_MOUSE
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081