

Human IL-12 receptor beta 1 protein (Recombinant) (STJP000542)

STJP000542

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human IL-12 receptor beta 1-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human Interleukin-12 receptor beta was lyophilized from a.2 Mu m filtered PBS solution pH7.4.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	3594
Gene Symbol	IL12RB1
Uniprot ID	I12R1_HUMAN
Immunogen Sequence	