

Human LIFR (CD118) protein (Recombinant) (STJP000543)

STJP000543

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human LIFR (CD118)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CD118 is supplied as a 0.2 µg/ml filtered PBS solution, pH7.2.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD118, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	3977
Gene Symbol	LIFR
Uniprot ID	LIFR_HUMAN
Immunogen Sequence	