

Mouse LIFR (CD118) protein (Recombinant) (STJP000544)

STJP000544

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse LIFR (CD118)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse CD118 (LIFR) is supplied as a 0.2 µg/ml filtered PBS solution, pH7.2.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC
Storage	Recombinant mouse CD118, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3977
Gene Symbol	LIFR
Uniprot ID	LIFR_HUMAN
Immunogen	
Sequence	