

Rat IL-1 receptor type II (CD121b) protein (Recombinant) (STJP000546)

STJP000546

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Rat IL-1 receptor type II (CD121b)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Rat IL-1 receptor II (CD121b) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage Instruction	Recombinant Rat IL1R2 can be stored in working aliquots at 2-8°C for one month, or at-20°C to-70°C for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	117022
Gene Symbol	Il1r2
Uniprot ID	IL1R2_RAT
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081