

## Mouse IL-12 protein (Recombinant) (STJP000547)

STJP000547

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse IL-12-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant mouse Interleukin-12 was lyophilized from a.2 Mu m filtered PBS solution pH7.4.
<b>Purification</b>	
<b>Dilution Range</b>	>90%, as determined by SDS-PAGE
<b>Storage</b>	The lyophilized protein is stable for at least years from date of receipt at -20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	16159
<b>Gene Symbol</b>	IL12a
<b>Uniprot ID</b>	IL12A_MOUSE
<b>Immunogen Sequence</b>	