

Human HSP47 protein (Recombinant) (STJP000558)

STJP000558

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human HSP47-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Heat shock protein 47 (HSP47) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Human heat shock protein 47 (SERPINH1) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	871
Gene Symbol	SERPINH1
Uniprot ID	SERPH_HUMAN
Immunogen Sequence	