

## Mouse CD38 protein (Recombinant) (STJP000562)

STJP000562

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse CD38-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant mouse CD38 is supplied as a 0.2 µg/ml filtered PBS solution, pH7.2.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant mouse CD38, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	12494
<b>Gene Symbol</b>	Cd38
<b>Uniprot ID</b>	CD38_MOUSE
<b>Immunogen</b>	
<b>Sequence</b>	