

Human BACE1 protein (Recombinant) (STJP000564)

STJP000564

GENERAL INFORMATION

Product Type Proteins
Short Description Recombinant-Human BACE1-protein was developed from hek293. For use in research applications.
Host/Source HEK293

PRODUCT PROPERTIES

Concentration
Formulation Recombinant Human BACE1 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification
Dilution Range >95%, as determined by SDS-PAGE and HPLC
Storage Recombinant BACE1 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid
Instruction repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID [23621](#)
Gene Symbol [BACE1](#)
Uniprot ID [BACE1_HUMAN](#)
Immunogen
Sequence