

Mouse CD80 (B7-1) protein (Recombinant) (STJP000583)

STJP000583

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse CD80 (B7-1)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse CD80 was lyophilized from 0.2um filtered 50mM Tris, 100mM NaCl, pH 8.0.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC
Storage	Recombinant mouse CD80, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	12519
Gene Symbol	Cd80
Uniprot ID	CD80_MOUSE
Immunogen	
Sequence	