

## Mouse PD-1 (PDCD1) protein (Recombinant) (STJP000588)

STJP000588

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse PD-1 (PDCD1)-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant Mouse PD-1 is lyophilized from 0.2 M filtered PBS solution, pH7.2, 5% Trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant PD-1 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	18566
<b>Gene Symbol</b>	Pdcd1
<b>Uniprot ID</b>	PDCD1_MOUSE
<b>Immunogen</b>	
<b>Sequence</b>	