

Human Glial growth factor (GGF2) protein (Recombinant) (STJP000590)

STJP000590

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Glial growth factor (GGF2)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant GGF2 was lyophilized from a.2 Mu m filtered0.02M PB, pH7.0.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	The lyophilized protein is stable for at least years from date of receipt at-20°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at-20°C for six months, with a carrier protein without detectable loss
Instruction	

TARGET INFORMATION

Gene ID	3084
Gene Symbol	NRG1
Uniprot ID	NRG1_HUMAN
Immunogen	
Sequence	