

## Rat IL-1 receptor 1 (CD121a) protein (Recombinant) (STJP000591)

STJP000591

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Rat IL-1 receptor 1 (CD121a)-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant rat IL-1R (CD121a) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>97%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant rat IL-1 receptor type 1 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">NA</a>
<b>Gene Symbol</b>	<a href="#">IL1r1</a>
<b>Uniprot ID</b>	<a href="#">IL1R1_RAT</a>
<b>Immunogen</b>	
<b>Sequence</b>	