

Human CD150 (SLAM) protein (Recombinant) (STJP000602)

STJP000602

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD150 (SLAM)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CD150 is supplied as a 0.2 µg/ml filtered PBS solution, pH7.2.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant SLAM (CD150), as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6504
Gene Symbol	SLAMF1
Uniprot ID	SLAF1_HUMAN
Immunogen	
Sequence	