

## Human FGFBP1 protein (Recombinant) (STJP000608)

STJP000608

### GENERAL INFORMATION

**Product Type** Proteins  
**Short Description** Recombinant-Human FGFBP1-protein was developed from hek293. For use in research applications.  
**Host/Source** HEK293

### PRODUCT PROPERTIES

**Concentration**  
**Formulation** Recombinant human Fibroblast growth factor binding protein (FGF-BP1) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.  
**Purification**  
**Dilution Range** >95%, as determined by SDS-PAGE and HPLC  
**Storage** Recombinant Human FGFBP1 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.  
**Instruction** Avoid repeated freeze/thaw cycles.

### TARGET INFORMATION

**Gene ID** [9982](#)  
**Gene Symbol** [FGFBP1](#)  
**Uniprot ID** [FGFBP1\\_HUMAN](#)  
**Immunogen**  
**Sequence**