

Human PD-1 (PDCD1) protein (Recombinant) (STJP000611)

STJP000611

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human PD-1 (PDCD1)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human PD-1 is lyophilized from 0.2 M unfiltered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant PD-1 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. It is
Instruction	recommended to add 0.1% HSA or BSA for long term storage. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	5133
Gene Symbol	PDCD1
Uniprot ID	PDCD1_HUMAN
Immunogen Sequence	