

## Human Adiponectin protein (Recombinant) (STJP000613)

STJP000613

### GENERAL INFORMATION


<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human Adiponectin-protein was developed from cho cells. For use in research applications.
<b>Host/Source</b>	CHO cells

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant Human Adiponectin is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>97%, as determined by SDS-PAGE and HPLC
<b>Storage Instruction</b>	Recombinant Adiponectin can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	9370
<b>Gene Symbol</b>	ADIPOQ
<b>Uniprot ID</b>	ADIPO_HUMAN
<b>Immunogen Sequence</b>	



Cells were treated with gAcrp (1 Mu g/ml) for the indicated time points.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081