

Human R-Spondin1 protein (Recombinant) (STJP000618)

STJP000618

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human R-Spondin1-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human R-spondin1 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage Instruction	Recombinant R-spondin1 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	284654
Gene Symbol	RSPO1
Uniprot ID	RSPO1_HUMAN
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081