

Rat CD155 protein (Recombinant) (STJP000623)

STJP000623

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Rat CD155-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant rat CD155is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant rat PVR (CD155) , as supplied, can be stored in working aliquots at 2-8°C for one month, or at-20°C to-70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	