

Mouse IL-18 receptor 1 protein (Recombinant) (STJP000627)

STJP000627

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse IL-18 receptor 1-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Mouse Interleukin-18 receptor alpha is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant IL-18 receptor alpha, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Sequence