

Mouse EPCR protein (Recombinant) (STJP000638)

STJP000638

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse EPCR-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Mouse Endothelial protein C receptor (PROCR) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Mouse EPCR can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	19124
Gene Symbol	Procr
Uniprot ID	EPCR_MOUSE
Immunogen Sequence	