

Human IL-20R2 protein (Recombinant) (Fc-Tag) (STJP000642)

STJP000642

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human IL-20R2-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human Interleukin-20 receptor beta (IL-20R2) is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant human IL-20RB, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	53833
Gene Symbol	IL20RB
Uniprot ID	I20RB_HUMAN
Immunogen Sequence	