

## Mouse PCSK9 protein (Recombinant) (STJP000647)

STJP000647

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse PCSK9-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant Mouse PSCK9 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>95%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant PCSK9 can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">100102</a>
<b>Gene Symbol</b>	<a href="#">Pcsk9</a>
<b>Uniprot ID</b>	<a href="#">PCSK9_MOUSE</a>
<b>Immunogen</b>	
<b>Sequence</b>	