

Mouse CD226 (DNAM-1) protein (Recombinant) (STJP000653)

STJP000653

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse CD226 (DNAM-1)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant mouse CD226 is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>97%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Mouse DNAM1 (CD226) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	225825
Gene Symbol	Cd226
Uniprot ID	CD226_MOUSE
Immunogen	
Sequence	