

Human Neuroligin-1 protein (Recombinant) (STJP000655)

STJP000655

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Neuroligin-1-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human Neuroligin-1 (NLGN1) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Neuroligin-1 (NLG1) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	57502
Gene Symbol	NLGN4X
Uniprot ID	NLGNX_HUMAN
Immunogen Sequence	