

Human IL-17RC protein (Recombinant) (Fc-Tag) (STJP000656)

STJP000656

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human IL-17RC-Fc-Tag protein was developed from hek293 and has a target region of Fc-Tag. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Interleukin-17RC is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant IL-17RC, as supplied, can be stored in working aliquots at 2-8°C for one month, or at-20°C to-70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	84818
Gene Symbol	IL17RC
Uniprot ID	I17RC_HUMAN
Immunogen Sequence	