

Human Scavenger Receptor Class B Type I protein (Recombinant) (STJP000663)

STJP000663

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Scavenger Receptor Class B Type I-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human SR-BI (SCARB1) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>90%, as determined by SDS-PAGE and HPLC
Storage	Recombinant SR-BI can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	949
Gene Symbol	SCARB1
Uniprot ID	SCRB1_HUMAN
Immunogen	
Sequence	