

## Human CD112 (PRR2) protein (Recombinant) (STJP000667)

STJP000667

### GENERAL INFORMATION

**Product Type** Proteins  
**Short Description** Recombinant-Human CD112 (PRR2)-protein was developed from hek293. For use in research applications.  
**Host/Source** HEK293

### PRODUCT PROPERTIES

**Concentration**  
**Formulation** Recombinant NECTIN2 was lyophilized from 0.2 µm filtered PBS solution pH7.4  
**Purification**  
**Dilution Range** >95%, as determined by SDS-PAGE and HPLC  
**Storage Instruction** The lyophilized protein is stable for at least 2 years from date of receipt at -20°C.

### TARGET INFORMATION

**Gene ID** [5819](#)  
**Gene Symbol** [NECTIN2](#)  
**Uniprot ID** [NECT2\\_HUMAN](#)  
**Immunogen Sequence**