

Human IL-17RA (CD217) protein (Recombinant) (STJP000671)

STJP000671

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human IL-17RA (CD217)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Interleukin-17 receptor alpha is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>98%, as determined by SDS-PAGE and HPLC
Storage	Recombinant IL-17 receptor alpha, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	23765
Gene Symbol	IL17RA
Uniprot ID	I17RA_HUMAN
Immunogen	
Sequence	