

Human CD114 (GCSFR) protein (Recombinant) (STJP000675)

STJP000675

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD114 (GCSFR)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant GCSFR (CD114) is supplied as a 0.2 Mu m filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD114 as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	1441
Gene Symbol	CSF3R
Uniprot ID	CSF3R_HUMAN
Immunogen Sequence	