

## Mouse Neuroligin-1 protein (Recombinant) (STJP000678)

STJP000678

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Mouse Neuroligin-1-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant Mouse Neuroligin-1 (NLGN1) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
<b>Purification</b>	
<b>Dilution Range</b>	>90%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant Neuroligin-1 (NLG1) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
<b>Instruction</b>	Avoid repeated freeze/thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">192167</a>
<b>Gene Symbol</b>	<a href="#">Nlgn1</a>
<b>Uniprot ID</b>	<a href="#">NLGN1_MOUSE</a>
<b>Immunogen Sequence</b>	