

Human Sclerostin protein (Recombinant) (STJP000679)

STJP000679

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | Recombinant-Human Sclerostin-protein was developed from hek293. For use in research applications. |
| Host/Source | HEK293 |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Concentration | |
| Formulation | Recombinant Human Sclerostin is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose. |
| Purification | |
| Dilution Range | >95%, as determined by SDS-PAGE and HPLC |
| Storage Instruction | Recombinant Sclerostin can be stored in working aliquots at 2-8°C for one month, or at-20°C to-70°C for twelve months. Avoid repeated freeze/thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|------------|
| Gene ID | 50964 |
| Gene Symbol | SOST |
| Uniprot ID | SOST_HUMAN |
| Immunogen Sequence | |