

Human CD113 (Nectin-3) protein (Recombinant) (STJP000694)

STJP000694

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human CD113 (Nectin-3)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant CD113 is supplied as a 0.2 µg/ml filtered PBS solution, pH7.2.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant CD113, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	25945
Gene Symbol	NECTIN3
Uniprot ID	NECT3_HUMAN
Immunogen Sequence	