

## Human IL-20RA protein (Recombinant) (STJP000704)

STJP000704

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Recombinant-Human IL-20RA-protein was developed from hek293. For use in research applications.
<b>Host/Source</b>	HEK293

### PRODUCT PROPERTIES

<b>Concentration</b>	
<b>Formulation</b>	Recombinant Interleukin-20 receptor alpha is supplied as lyophilized 0.2 Mu m filtered PBS solution, pH7.2.
<b>Purification</b>	
<b>Dilution Range</b>	>98%, as determined by SDS-PAGE and HPLC
<b>Storage</b>	Recombinant IL-20 receptor alpha, as supplied, can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	53832
<b>Gene Symbol</b>	IL20RA
<b>Uniprot ID</b>	I20RA_HUMAN
<b>Immunogen</b>	
<b>Sequence</b>	