

Human Kallikrein-4 protein (Recombinant) (STJP000711)

STJP000711

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Kallikrein-4-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Kallikrein-4 (KLK-4) is lyophilized from 0.2 Mu m filtered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant Kallikrein-4 (KLK4) can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months.
Instruction	Avoid repeated freeze/thaw cycles.

TARGET INFORMATION

Gene ID	9622
Gene Symbol	KLK4
Uniprot ID	KLK4_HUMAN
Immunogen Sequence	