

Mouse Interleukin-18 Binding Protein (IL-18BP) protein (Recombinant) (STJP000712)

STJP000712

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Mouse Interleukin-18 Binding Protein (IL-18BP)-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant IL-18BP is supplied as an.2 Mu m filtered PBS solution.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Upon reconstitution, this cytokine can be stored in working aliquots at 2-8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	16068
Gene Symbol	Il18bp
Uniprot ID	I18BP_MOUSE
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081