

Human EGFR (erbb1) protein (Recombinant) {Biotin} (STJP000714)

STJP000714

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human EGFR (erbb1)-protein was developed from cho cells. For use in research applications.
Host/Source	CHO cells

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human EGFR was lyophilized from a 0.2 M filtered PBS pH 7.4 with 10% trehalose
Purification	
Dilution Range	>95%, as determined by SDS-PAGE
Storage	Recombinant EGFR can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. It is recommended to add 0.1% HSA or BSA for long term storage. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1956
Gene Symbol	EGFR
Uniprot ID	EGFR_HUMAN
Immunogen Sequence	