

Human Lipocalin-2 protein (Recombinant) (STJP000744)

STJP000744

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant-Human Lipocalin-2-protein was developed from hek293. For use in research applications.
Host/Source	HEK293

PRODUCT PROPERTIES

Concentration	
Formulation	Recombinant Human Lipocalin-2 is lyophilized from 0.2 M unfiltered PBS solution, pH7.2, 5% Trehalose.
Purification	
Dilution Range	>95%, as determined by SDS-PAGE and HPLC
Storage	Recombinant NGAL can be stored in working aliquots at 2-8°C for one month, or at -20°C to -70°C for twelve months. Avoid repeated freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3934
Gene Symbol	LCN2
Uniprot ID	NGAL_HUMAN
Immunogen	
Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081