

Mouse Anti-HBV HBcAg antibody {Gold} (STJS000084)

STJS000084

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Gold-conjugated mouse monoclonal anti-HBV HBcAg secondary antibody. For use in most research applications.
Applications	WB, IHC-P, IF, ICC, ELISA
Host/Source	Mouse
Reactivity	HBV

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	4-5 mg/mL
Conjugation	Gold
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:5000-10000 IF 1:5000-10000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence