

Mouse Anti-Pig IgG Heavy & Light Chain antibody {HRP} (STJS000121)

STJS000121

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated mouse monoclonal anti-Pig IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	ELISA, WB, IHC, ICC
Host/Source	Mouse
Reactivity	Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	HRP
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	ELISA 1:5000-1:20000 WB 1:5000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	