

Rabbit Anti-Rat IgG Heavy & Light Chain antibody {HRP} (STJS000138)

STJS000138

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated rabbit polyclonal anti-Rat IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	ELISA, WB, IHC, ICC
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	HRP
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	ELISA 1:5000-1:20000 WB 1:5000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence