

Rabbit Anti-Rat IgG Heavy & Light Chain antibody {Biotin} (STJS000170)

STJS000170

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Secondary antibodies |
| Short Description | Biotin-conjugated rabbit polyclonal anti-Rat IgG Heavy & Light Chain secondary antibody. For use in most research applications. |
| Applications | ELISA, WB, IHC, ICC |
| Host/Source | Rabbit |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Biotin |
| Purification | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Dilution Range | ELISA 1:5000-1:20000 WB 1:5000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000 |
| Formulation | 0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.05% Sodium Azide |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |