

Mouse Anti-Monkey IgG Heavy & Light Chain antibody {Biotin} (STJS000174)

STJS000174

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Biotin-conjugated mouse monoclonal anti-Monkey IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	ELISA, WB, IHC, ICC
Host/Source	Mouse
Reactivity	Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Biotin
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	ELISA 1:10000-1:20000 WB 1:10000-1:20000 IHC 1:500-1:1000 ICC 1:500-1:1000
Formulation	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence