

Rabbit Anti-Guinea Pig IgG Heavy & Light Chain antibody {FITC} (STJS000193)

STJS000193

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	FITC-conjugated rabbit polyclonal anti-Guinea Pig IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	IF, ICC, FC
Host/Source	Rabbit
Reactivity	Guinea Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	FITC
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IF 1:250-1:2000 FCM 1:250-1:2000
Formulation	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	