

**Mouse Anti-Duck IgY Heavy & Light Chain antibody {FITC} (STJS000205)**

STJS000205

**GENERAL INFORMATION**

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	FITC-conjugated mouse monoclonal anti-Duck IgY Heavy & Light Chain secondary antibody. For use in most research applications.
<b>Applications</b>	IF, ICC, FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Duck

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	FITC
<b>Purification</b>	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	IF 1:250-1:2000 FCM 1:250-1:2000
<b>Formulation</b>	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.05% Sodium Azide
<b>Isotype</b>	IgY
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>