

## Rabbit Anti-Goat IgG Heavy & Light Chain antibody {Alexa Fluor 350} (STJS000225)

STJS000225

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Secondary antibodies  |
| <b>Short Description</b> | Alexa Fluor 350-conjugated rabbit polyclonal anti-Goat IgG Heavy & Light Chain secondary antibody. For use in most research applications. |
| <b>Applications</b>      | ELISA, IF, ICC, FC  |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Goat  |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 1 mg/mL   |
| <b>Conjugation</b>         | Alexa Fluor 350   |
| <b>Purification</b>        | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| <b>Dilution Range</b>      | IHC 1:200-1:1000<br>IF 1:200-1:1000<br>FCM 1:2000-1:4000<br>ELISA Use at an assay dependent concentration         |
| <b>Formulation</b>         | 0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                    |

### TARGET INFORMATION

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            |  |
| <b>Gene Symbol</b>        |  |
| <b>Uniprot ID</b>         |  |
| <b>Immunogen</b>          |  |
| <b>Immunogen Region</b>   |  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |