

**Mouse Anti-Human Ig Lambda Chain antibody {Alexa Fluor 350} (STJS000248)**

STJS000248

**GENERAL INFORMATION**

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	Alexa Fluor 350-conjugated mouse monoclonal anti-Human Ig Lambda Chain secondary antibody. For use in most research applications.
<b>Applications</b>	ELISA, IF, ICC, FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Alexa Fluor 350
<b>Purification</b>	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	IHC 1:200-1:1000
	IF 1:200-1:1000
	FCM 1:2000-1:4000
	ELISA Use at an assay dependent concentration
<b>Formulation</b>	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide
<b>Isotype</b>	Ig
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>