

## Mouse Anti-Human Ig Lambda Chain antibody {Alexa Fluor 405} (STJS000291)

STJS000291

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Secondary antibodies  |
| <b>Short Description</b> | Alexa Fluor 405-conjugated mouse monoclonal anti-Human Ig Lambda Chain secondary antibody. For use in most research applications. |
| <b>Applications</b>      | ELISA, IF, ICC, FC  |
| <b>Host/Source</b>       | Mouse   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            |  |
| <b>Concentration</b>       | 1 mg/mL  |
| <b>Conjugation</b>         | Alexa Fluor 405  |
| <b>Purification</b>        | The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads. |
| <b>Dilution Range</b>      | IHC 1:200-1:1000<br>IF 1:200-1:1000<br>FCM 1:2000-1:4000<br>ELISA Use at an assay dependent concentration              |
| <b>Formulation</b>         | 0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide  |
| <b>Isotype</b>             | Ig   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                         |

### TARGET INFORMATION

|                         |
|-------------------------|
| <b>Gene ID</b>          |
| <b>Gene Symbol</b>      |
| <b>Uniprot ID</b>       |
| <b>Immunogen</b>        |
| <b>Immunogen Region</b> |
| <b>Specificity</b>      |
| <b>Immunogen</b>        |
| <b>Sequence</b>         |